




## Product Data Sheet



Rev\_061316

# HCMV Grade 2 Antigen Solution (ready-to-use)

**Catalog number:** HCMV-GR2-015

**Product Name:** HCMV Grade 2 Antigen Solution, inactivated by gamma radiation

**Package Size:** HCMV Grade 2 antigen dissolved in serum-free CTL-Test™ Medium, 1.5ml per vial

**Concentration:** 60µg/ml

**Description:** This preparation contains HCMV antigens from all parts of the virus cycle, and residual host cellular material, (inactivated by gamma radiation).

**Activity:** HCMV stimulates antigen-specific CD4+ memory T cells.

**Recommended Use:** After thawing, the reagent is a 2x ready-to-use working solution (concentration 60 µg/ml). When working in 96 well plates, add 100µL of antigen solution to 100uL of PBMC per well. The PBMC concentration should be 100,000 – 800,000 cells/well.

**Storage:** Keep frozen at -80°C.

**For laboratory research use only. Not for use in diagnostic or therapeutic procedures.**